



# Sonic Denture Cleaner

## INSTRUCTION MANUAL

Model No. EL12378 ((XL-E0005)

Please read carefully before use and keep for future reference

A major advance in denture hygiene, this battery-powered cleaner uses sonic vibrations to help remove bacteria, stains, plaque and food particles. Similar to the type employed by dentists and dental hygienists, it's wonderfully easy to use. Just pop your dentures in the container along with your usual cleaning solution. Switch on and in only a few minutes you'll appreciate how much cleaner, whiter and fresher they are. And no brushing required!

### **Operation:**

Pull out the battery level at the bottom of the denture cleaner and insert 2 x AA batteries (not included) observing correct polarity as indicated in the battery compartment.

### Dentures

- Remove the container from your Sonic Cleaner and place your dentures inside.
- Cover the dentures with your standard cleaning solution. Always follow your dentures manufacturer instructions.  
**DO NOT USE BOILING WATER.**  
(Use with any denture, retainer cleaner or just fresh tap or tepid water)
- Turn the switch located at back side of the cleaner to ON position and let the machine work for 5 minutes.
- Please note you can run the machine for longer, but it is better to have a pause every 5 minutes to check the condition of your dentures.
- Switch off after 5 minutes, rinse the dentures in fresh tap or tepid water and brush lightly with a soft bristle brush.

### **Jewellery**

- Place the jewellery on the tray.
- Cover the jewellery with a suitable cleaning solution or some fresh tap or tepid water. Do not use any cleaning foaming agents.
- Turn the switch located at back side of the cleaner to ON position and let the machine work for 5 minutes. Repeat if necessary.

**THIS PRODUCT IS NOT RECOMMENDED FOR PEARLS, PEARL COATED JEWELLERY, ABALONE OR ANY TYPE MOLLUSC NACRE JEWELLERY.**

### **Technical Specification:**

Container holds 150cc.

Requires 2 x AA batteries (not included).

Measurement: 9 x 8 x 7.5cm

Imported and distributed by EasyLife Group Ltd, London NW1 0JP

**Support Line: 0303 031 0777**

For more innovative products from EasyLife please visit:

**[www.easylifegroup.com](http://www.easylifegroup.com)**